

LABORATORY SERVICES

TONHC Specimen Collection Overview - 02.18.2020







SPECIMEN LABELING



Proper labeling of specimens

FULL

PATIENT NAME, DOB, CHART NUMBER

FULL

PATIENT NAME, DOB

FULL

PATIENT NAME, CHART NUMBER


- Per NPSG 01.01.01:

 labels must have at least
 patient identifiers AND
 be labeled while patient
 is present
- Incorrectly labeled specimens will not be processed - please use caution
- Culture labels should include site of collection



Lab Specimens





- Clean catch urine
- Protein reflex
- Pregnancy test
- Urine GC/Chlamydia
- Urine HCG



Common Swabs



Rapid Strep – Double swab



- Influenza
 Nasopharyngeal Swab
 (A or B)
 - Send 2 swabs in any combination
 - Please note these are FOAM tipped only




- Group Beta Strep
- Wound Culture
- Fecal Leukocytes



Viral testing





- RSV (items given by lab)
- RED TOP TUBE & WIRE SWAB

Viral Media/HSV



Stool Specimens




- Ova & Parasite Exam
- Giardia



Stool culture



Additional Stool Tests



- Clostridium difficile (C. Diff)
- Helicobacter pylori (H. Pylori)



Additional Stool Tests



Hemoccult In-house test



Occult Blood-FIT
Send-out test



^{*} Each test has its own unique order in CPRS. Order entry must match slide collected.

Blue tube specifics



- BOTH tubes exist in the TONHC inventory
- Either tube can be used when a "BLUE TUBE" is indicated for a specimen
- Please fill to line indicated by black arrows.
- The tube with the lower fill line will also have a white ring on the top aspect



Blood Cultures



One pair of bottles per Blood Culture order if adult.
One pink bottle for Pediatric.



Adult specimen

Pediatric or younger





LAB FORMS COMPLETION



Newborn Screening

SPECIMEN	Date / Ti	me Stamp		00
aby's Name ast: First:			Submitter / Physician Information AZ252317581	5
Date of Birth	Time of Birth Birth	Weight Sex Grams M	Submitter Name/ID:	1 12 1
Date of Collection	Time of Collection Curre		Ordering Physician (Hosp.):Follow-up Physician Name (Last,First):	23 Ahls
Baby's AHCCCS #		Gestational Age	Phone: ()	5
MR #		WeeksDays	Practice Address:	NI
	Multiple Birth (circle one) A	B C D	City, State, Zip:	8 N
Race 1 White 2 African Amer.	Food Source 1 Breast Only 2 Milk / Lactose Formula	Status Meconium Ileus Y N In NICU/Special Y N	Birth Mother's Information Mom's Name Last:First:	2020-11 26 A
3 Asian 4 Amer. Indian 5 Other	3 Soy Formula 4 Breast & Lactose	Transfusion (RBC ONLY) before collection?	Mom's Date of Birth: / Maiden Name: (OR) Other Person with Custody:	er 2.
Hispanic Y N	6 TPN 0 Never Fed	Date FIRST transfused	Street Address: Insurance papers included	EXP DATE PerkinElme Lot 1056



SF-518

518-123				NSN 7540-00-634-4158						
MEDICAL RECORD		BLOOD OR	BLOOD COMPONENT	TRANSFUSION						
	•		REQUISITION							
COMPONENT REQUESTED	(Check one)	TYPE OF REQUEST (Chec. Products are requested.)	k ONLY if Red Blood Cell	REQUESTING PHYSICIAN (Print)						
RED BLOOD CELLS		TYPE AND SCREEN		DI ONOCIO OD ODERATIVE DROCEDINE						
FRESH FROZEN PLASMA				DIAGNOSIS OR OPERATIVE PROCEDURE						
PLATELETS (Pool of	units)	CROSSMATCH								
CRYOPRECIPITATE (Pool of units)		DATE REQUESTED		I have collected a blood specimen on the below named						
Rh IMMUNE GLOBULIN OTHER (Specify) VOLUME REQUESTED (if applicable)		DATE AND HOUR REQUIRED KNOWN ANTIBODY FORMATION/TRANSFUSION REACTION (Specify)		patient, verified the name and ID No. of the patient and verified the specimen tube label to be correct. SIGNATURE OF VERIFIER						
						ML		TABLE TO THE (OPERATO)		
						REMARKS:		IF PATIENT IS FEMALE, IS	THERE HISTORY OF:	DATE VERIFIED
		RhIG TREATMENT? DATE GIVEN:		TIME VERIFIED						
		HEMOLYTIC DISEASE OF NEWBORN?		THE VERNIED						
		SECTION II - PRE-TR.	ANSFUSION TESTING							
UNIT NO.	TRANSFUSION NO.		RPRETATION	PREVIOUS RECORD CHECK:						
	PATIENT NO.	ANTIBODY SCREEN	CROSSMATCH	RECORD NO RECORD SIGNATURE OR PERSON PERFORMING TEST						
	PATIENT NO.			SIGNATURE SITT ENGINE ENGINEERS						
DONOR	RECIPIENT									
ABO	ABO	REMARKS:	EQUIRED FOR THE COMPO	NENT REQUESTED DATE						
200	200	NEWSTAND.								
Rh	Rh									
		SECTION III DECOR	RD OF TRANSFUSION							
	PRE-TRANSFUSION DATA	SECTION III - RECOR	T TRANSPOSION	POST-TRANSFUSION DATA						
INSPECTED AND ISSUED B			AMOUNT GIVEN	TIME/DATE COMPLETED/INTERRUPTED						
AT (Hour)	ON (Date)		REACTION NONE SUSPECTED	TEMPERATURE PULSE BLOOD PRESSURE						
IDENTIFICATION	011 (2010)		If reaction is suspected – IMI							
I have examined the Blood information identifying the co			Discontinue transfusion, tr Notify Physician and Trans	reat shock if present, keep intravenous line open. sfusion Service.						
	rson named on this Blood C	component Transfusion Form	Follow Transfusion Reacti							
1st VERIFIER (Signature)			DESCRIPTION OF REACTION							
			URTICARIA CH	HILL FEVER PAIN						
2nd VERIFIER (Signature)			OTHER (Specify)							
			OTHER RIESION TIES (5							
PRE-TRANSFUSION			OTHER DIFFICULTIES (Equ.	upment, clots, etc.)						
TEMP.	PULSE	BP	SIGNATURE OF PERSON N	NOTING ABOVE						
DATE OF TRANSFUSION	TIME STARTED)								
PATIENT IDENTIFICATION -	 USE EMBOSSER (For typed rate; hospital or medical facil 	l or written entries give: Name– ity)	Last, first, middle; grade; rank;	SEX WARD						
				BLOOD OR BLOOD COMPONENT TRANSFUSION Medical Record						
				STANDARD FORM 518 (REV. 9-92) Prescribed by GSA/ICMR, FIRMR (41 CFR) 201-9.202-1						



Blood Consent

Tohono O'odham Nation Health Care Sells Hospital Transfusion of Blood/Blood Products Consent Form

Consen		
	blood/blood product(s nd his/her associates :	
(Specify type of product)		
2. The reason for the transfusion has been explained	I to me as	
(Layman terms) Alternatives to transfusions may not always work when much care is taken, transfusion can cause p		ition. I also understand that even
 a.) Allergic Reactions b.) Hepatitis c.) Acquired Immune Deficiency Syndrome (d.) Develop allergies that will make future trae.) Possible exposure to other diseases. 		cult.
 The above has been explained to me by my phys and all my questions have been answered. 	ician(Print phy	ysician name)
4. I have been informed of alternatives to receiving	a blood transfusion/b	lood products.
5. I accept on behalf of myself (the patient) the risk	s associated with a tra	unsfusion as described above.
Patient Signature		Date Date
Patient Representative/Interpreter (if patient is unable to sign)		Date
Physician Signature	Time	Date
Witness Signature		Date
ime:	1	
DB:		
nart#:		
	I	TONHC Pt Consent

Updated 1.2020 tma



Emergency

Tohono O'odham Nation Health Care Sells Hospital P.O. Box 548 Highway 86 Sells, Arizona 85634

******* <u>UNCROSSMATCHED</u> *******

EMERGENCY RELEASE OF BLOOD COMPONENTS WITHOUT COMPLETION OF BLOOD BANK TESTS

I, the undersigned, have deemed the clinical situation was sufficiently emergent to require release of uncrossmatched blood components. I affirm the benefit of using uncrossmatched blood outweighs the risk. I assume all responsibility for ordering blood (packed red blood cells) to be given to the patient, named below, without completion of blood bank compatibility testing. I further understand that the blood bank laboratory technical staff will perform compatibility testing as soon as possible and that any incompatible test result will be reported to me immediately. In my professional medical judgment, this patient's life will be in jeopardy without an emergency transfusion. units of: Group O, Rh negative red blood cells ☐ ABO type specific red blood cells ☐ Other, specify: ____ or, R.N.: M.D. Date: Time: For M.D.: Lab Use Only: Tech Initials ____ PATIENT IDENTIFICATION - USE EMBOSSER

> Medical Records TONHC - 1BB Implemented 4/2017

Document or adhere BUI sticker(s) of units here.






CRITICAL LAB REPORTING



Critical Lab Reporting


- NPSG.02.03.01 requires reporting of critical diagnostic test results on a timely basis
- Per TONHC policy it is mandatory to complete this within 15 minutes.





Recording Critical Labs

- Activate each section by clicking "Critical Lab Values #(X)"
- Complete one section for each lab value





Attention to detail can't be (and never is) added later.

- Marco Arment



QUESTIONS?





TONHC Sells Laboratory Services

520-383-7233

